

BIOMECHANICAL STUDY SUMMARY

Indirect Decompression and Vertebral Body Endplate Strength after Lateral Interbody Spacer Impaction: Cadaveric and Foam-Block Models

Anthony J. Kwon, William D. Hunter, Mark Moldavsky, Kanaan Salloum, and Brandon Bucklen

Journal of Spine Surgery 4(1):62-71, 2018.


OBJECTIVE: The primary aim of this study was to determine whether implant trialing and impaction of a static lateral spacer results in reduced endplate strength compared to that achieved by insertion of a CALIBER®-L expandable lateral spacer.

METHOD: Using a cadaveric biomechanical and foam-block vertebral model, researchers compared vertebral body endplate strength and distraction potential between static and CALIBER®-L expandable spacer groups. Fourteen lumbar motion segments (seven L2-3 and seven L4-5 specimens) were distributed evenly between static and CALIBER®-L expandable spacer groups. In each specimen, discectomy was followed by trialing and spacer impaction. Motion segments were axially sectioned through the disc, and a metal stamp was used to apply a compressive load to superior and inferior vertebral bodies to quantify endplate strength. A paired, 2-sample for means t-test was performed to determine statistically significant differences between groups ($P \leq 0.05$). A foam-block endplate model was used to control simulated disc tension when a spacer with 2mm and 3mm desired distraction was inserted. One-way ANOVA and a *post hoc* Student Newman-Keuls test were performed ($P \leq 0.05$) to determine differences in distraction.



RESULTS:

- Both static and CALIBER®-L expandable spacers restored neural foraminal and disc height to intact levels, after discectomy and device implantation ($P > 0.05$)
- Maximum peak loads at endplate failure for static and CALIBER®-L expandable spacers were 1764N (± 966 N) and 2284N (± 949 N), respectively ($P \leq 0.05$).
- CALIBER®-L expandable spacers consistently produced greater desired distraction than that created by the static spacers in the foam-block model ($P \leq 0.05$).
- Distraction created by fully expanded CALIBER®-L expandable spacers was significantly greater than the predetermined goals of 2mm and 3mm ($P \leq 0.05$).



Measured versus desired distraction of static and CALIBER®-L expandable spacers at 240N resistance. * $P \leq 0.05$.

CONCLUSION: This biomechanical study showed that the increased trialing required for a static spacer may lead to additional iatrogenic endplate damage, resulting in less distraction and increased propensity for postoperative implant subsidence secondary to endplate disruption.



Scan the QR code to download the article or visit globusmedical.com/resource/caliber-l-kwon-et-al-2016

Talk to your Globus Medical sales representative to learn more about our complete line of expandable devices.

GMSS72
11.20 Rev A



GLOBUS
MEDICAL