DI110A
(REV. B)

ALLIANCE™ POSTERIOR LUMBAR DISC

GLOBUS MEDICAL, INC.
Vertical Force Business Center
2560 General Armistead Avenue
Audubon, PA 19403
USA
Customer Service:
Phone 1-866-GLOBUS1 (OR)
1-866-456-2871
Fax 1-866-GLOBUS3 (OR)
1-866-456-2873

IMPORTANT INFORMATION ON THE ALLIANCE™ POSTERIOR LUMBAR DISC
WICHTIGE INFORMATION ZUR POSTERIOREN LENDBANDSCHIEBE ALLIANCE™
INFORMATIONS IMPORTANTES SUR LE DISQUE LOMBAIRE POSTERIEUR ALLIANCE™
INFORMAZIONI IMPORTANTI SUL DISCO POSTERIORE LOMBARA ALLIANCE™
INFORMACIÓN IMPORTANTE SOBRE EL DISCO LUMBAR POSTERIOR ALLIANCE™
INFORMAÇÕES IMPORTANTES SOBRE O ALLIANCE™ POSTERIOR LUMBAR DISC
BELANGRIKKE INFORMATIE OVER DE ALLIANCE™ POSTERIOR LUMBARA DISCUS
VITKIĆ INFORMACIJE OZIJE ALLIANCE™ POSTERIOR LUMBAR DISK
ΤΑΡΚΕΑ ΤΕΤΟΑ ALLIANCE™ POSTERIORINEN LÄNNERANGAN VALIΛΕΥΡΠΡΟΤΕΣΙ
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ ΓΙΑ ΤΟΝ ΟΠΙΩΘΙΟ ΟΣΦΥΓΚΟ ΔΙΣΚΟ ALLIANCE™
ALLIANCE™ POSTERIOR LUMBAR DISC についての重要な情報

[EU REP] RMS - UK Limited
28 Trinity Road, Nailsea, Somerset, BS484NU England

CE 0297

FOR OUTSIDE THE UNITED STATES ONLY

ENGLISH

IMPORTANT INFORMATION
ON THE ALLIANCE™ POSTERIOR LUMBAR DISC

DESCRIPTION

The ALLIANCE™ Posterior Lumbar Disc is an implant assembly inserted using a posterior lumbar approach. The assembly consists of two implant assemblies, each with a superior and inferior endplate with keels and a moving polyethylene core. Implants are offered in a variety of configurations to accommodate varied patient anatomy. The ALLIANCE™ Posterior Lumbar Disc provides 0°, 4°, 8° or 12° of motion. Disc heights range from 1 mm to 16 mm. Two implant footprints are available with either 22mm or 25mm total implanted width.

The ALLIANCE™ endplates are manufactured from cobalt-chrome-molybdenum alloy, CoCrMo, as specified in ASTM F1537 (and ISO 5832-12). The superior and inferior surfaces of the ALLIANCE™ endplates are plasma sprayed with commercially pure titanium, as specified in ASTM F1580, F1978, F1147 and C-633 (and ISO 5832-2).

The ALLIANCE™ cores are manufactured from ultra-high molecular weight polyethylene, UHMWPE, as specified in ASTM F648 (and ISO 5834-2).

INDICATIONS

The ALLIANCE™ Posterior Lumbar Disc is indicated for spinal arthroplasty in treating skeletally mature patients with degenerative disc disease, primary or recurrent disc herniation, spinal stenosis, or spondylosis in the lumbosacral spine (L1-S1). DDD is defined as discogenic back pain with degeneration of the disc confirmed by patient history and radiographic studies, with or without leg (radicular) pain. Patients may have spondylolisthesis up to Grade 1 at the involved level.

CONTRAINdications

The ALLIANCE™ Posterior Lumbar Disc is contraindicated for patients with the following:

- Prior fusion surgery adjacent to the vertebral level being treated
- Prior surgery at the level to be treated
- Clinically compromised vertebral bodies at the affected level(s) due to current or past trauma
- Radiographic confirmation of facet joint disease or degeneration
- Osteoporosis, osteopenia, Paget's disease, osteomalacia or any other metabolic bone disease
- Active systemic or local infection
- Demonstrated allergy or foreign body sensitivity to any of the implant materials
- Conditions that may place excessive stress on the lumbar spine and implant, such as severe obesity or degenerative disease
- Patients whose activity, mental capacity, mental illness, alcoholism, drug abuse, occupation, or lifestyle may interfere with their ability to follow postoperative restrictions and who may place undue stresses on the implant during healing and may be at a higher risk of implant failure

WARNINGS

Risks associated with this device and with the surgical procedure are not known but may include:

- Injury to nerves, vessels, and organs
- Allergic or tissue reaction to implant materials
- Hematoma and impaired wound healing
- Migration or loosening of the device
- Fracture of the device
- Heterotopic ossification and fusion
- Venous thrombosis, lung embolism, and cardiac arrest
- Death

PRECAUTIONS

The implantation of spinal devices should be performed only by experienced spinal surgeons with specific training in the use of this system because this is a technically demanding procedure presenting a risk of serious injury to the patient. Preoperative planning and patient anatomy should be considered when selecting implant size.

ATTENTION

Under no circumstances should the ALLIANCE™ device components be used in combination with components from other suppliers. Damage to weight bearing structures can give rise to loosening of the device components, dislocation, and migration, and other serious complications; the implant should be checked periodically.

STERILIZATION

Implants:

ALLIANCE™ implants are supplied STERILE. The implant components may not be resterilized for any reason.

Instruments:

These instruments are supplied NONSTERILE. Sterilization is recommended as follows:

Method	Cycle	Temperature	Exposure Time
Steam	Gravity Displacement (Wrapped)	132° - 135°C (270° - 275°F)	20 Minutes
Steam	Pre-vacuum (Wrapped) Preconditioning Pulses: 3	132° - 135°C (270° - 275°F)	4 Minutes

These parameters are validated to sterilize only these instruments. If other products are added to the sterilizer, the recommended parameters are not valid and new cycle parameters must be established by the user. The autoclave must be properly installed, maintained, and calibrated. Ongoing testing must be performed to confirm inactivation of all forms of viable microorganisms.

SYMBOL TRANSLATION	
QTY	QUANTITY

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