<b>DI171B (REV. B)</b>	<b>KINEX® BIOACTIVE</b>
6/2024	
	<b>IMPORTANT INFORMATION ON KINEX® BIOACTIVE</b> WICHTIGE INFORMATIONEN ZU KINEX® BIOACTIVE INFORMATIONS IMPORTANTES CONCERNANT LE PRODUIT KINEX® BIOACTIVE INFORMAZIONI IMPORTANTI SU KINEX® BIOACTIVE INFORMACIÓN IMPORTANTE SOBRE KINEX® BIOACTIVE INFORMAZIONE IMPORTANTI ACRICA KINEX® BIOACTIVE WIKIG INFORMATION ON KINEX® BIOACTIVE WICHTIGE OPLYSNINGER OM KINEX® BIOACTIVE INFORMAZIONE IMPORTANTI OVER KINEX® BIOACTIVE ЗНАКОВАТА ИНФОРМАЦИЯ ЗА ПРОДУКТА KINEX® BIOACTIVE СЪЩЕСТВЕНА ИНФОРМАЦИЯ ЗА ПРОДУКТА KINEX® BIOACTIVE
<b>GLOBUS MEDICAL, INC.</b> Valley Forge Business Center 2500 General Armstrong Avenue Aurizton, PA 19403	 0297
Customer Service Phone 1-866-GLOBUS1 (CR) 1-866-456-2871 Fax 1-866-GLOBUS3 (CF) 1-866-456-2873	

KINEX® Bioactive is a resorbable bone void filler for the repair of bony defects. It is an osteoconductive and osteostimulative material that guides bone regeneration. When KINEX® is placed in direct contact with host bone, new bone grows in apposition to the surfaces of the implant. As the implant resorbs, bone and other connective tissues grow into the space previously occupied by KINEX®.

KINEX® implants consist of bings (per ASTM F1538), Type I bovine collagen (per ASTM F2212), and hyaluronic acid, and are available in gel, and strip forms to accommodate surgical and anatomical needs.

**INDICATIONS**  
KINEX® Bioactive is intended for use as a bone void filler and autograft extender for voids or gaps that are not intrinsic to the stability of the bony structure. These osseous defects may be surgically created or caused from traumatic injury to the bone. KINEX® is intended to be gently packed into bony voids or gaps of the skeletal system i.e., the extremities, pelvis, and spine) and should be combined with bone marrow aspirate. KINEX® resorbs and is replaced with bone during the healing process.

**WARNINGS**  
KINEX® is not designed with sufficient mechanical strength to support reduction of a defect site prior to soft and hard tissue ingrowth. Rigid fixation methods are recommended as needed to ensure stabilization of the defect. Complete postoperative wound closure is essential.

KINEX® Putty is intended for manual application and is not intended for injection through a constrained opening or under high pressure. KINEX® Gel may be injected into the desired location. High pressure injection of KINEX® should not be conducted as pressurization in closed cavities can lead to device extrusion beyond the intended application site, which could damage surrounding tissues, or to embolization of fat to the device into the bloodstream.

Packaging should be intact upon receipt. Damaged packages and/or contaminated products should not be used and should be returned to Globus Medical.

**PRECAUTIONS**  
KINEX® Bioactive is intended for use by surgeons familiar with bone grafting techniques. If fixation is used, the labeling for the use of the fixation system chosen should be followed and the fixation must gain purchase in the host bone. Standard postoperative practices for the treatment and rehabilitation associated with bone grafting must be strictly followed.

Mental or physical impairment which compromises a patient's ability to comply with necessary limitations or precautions may place a patient at a particular risk during postoperative rehabilitation.

A successful result is not always achieved in every surgical case. This is particularly true in spinal surgery where many extrinating circumstances may compromise the results.

As with any postoperative, care should be demonstrated in treating patients with preexisting conditions that may impact the success of the surgical procedure. This includes patients with bleeding disorders of any type, long-term steroid therapy, or immunosuppressive therapy, or high dose radiation therapy.

Use this device as supplied and in accordance with the handling and use information provided below.

**CONTRAINDICATIONS**  
KINEX® must not be used in patients with a history of anaphylaxis, history of multiple allergies, known allergies to bovine collagen who are being treated for desensitization to meat products because this product contains bovine collagen.

Conditions representing relative contraindications include:

- Severe neurological or vascular disease
- Hypocalcemia
- Pregnancy
- Cases of fracture fixation or where load support is required, unless standard internal or external stabilization techniques are followed to obtain rigid stabilization in all planes.
- Systemic and/or metabolic disorders that affect the bone or wound healing
- Conditions in which general bone grafting is not advisable
- Local infection
- Any patient unwilling to follow postoperative instructions
- Vertebral injury or kyphoplasty procedures
- To gain screw purchase or to stabilize screw placement
- Any case not described in the indications

#### POTENTIAL ADVERSE EVENTS

Possible complications are the same as to be expected of autogenous bone grafting procedures and include but are not limited to:

- Deformity of the bone at the surgical site
- Fracture or extrusion of the KINEX® implant(s), with or without generation of particulate debris
- Wound complications including hematoma, site drainage, infection (superficial, deep or deep with osteomyelitis), bone fracture, and other complications common to any surgical procedure
- Incomplete, or lack of, osseous ingrowth into bone void, as possible with any bone void filler
- Delayed union or failure of fusion
- Transient hypercalcemia
- Loss of bone graft, graft protrusion and/or dislodgement, which could damage surrounding tissue
- General complications that may arise from anesthesia and/or surgery

Localized immunological reactions consisting of transient localized edema, swelling, and rash have been reported to occur with bone void fillers containing collagen. Although there is evidence that the device will be assimilated or reabsorbed in some patients, the safety and effectiveness of the device in these patients has not been established.

Occurrence of one or more of these conditions may require an additional surgical procedure and may also require removal of the bone void filler.

**OSTEOSTIMULATION**  
KINEX® is an osteostimulative and osteoconductive device. Osteostimulation is defined as an accelerated bone formation process, characterized by the active stimulation of osteoblast proliferation and differentiation in an osseous defect. This stimulatory action has been demonstrated during *in vivo* tests of Bioglass particles to be more rapid than simple osteoconduction.<sup>1,2</sup> These tests have been supported by *in vivo* cell culture tests, which demonstrate the mechanisms of cell migration as being the result of cellular interaction with the ionic dissolution products released from Bioglass particles during its absorption.<sup>3-4</sup> Clinical data on this acceleration of bone formation in the human has not been established.

KINEX® has only been demonstrated to form bone in osseous defects. KINEX® therefore is not osteoconductive. Such osteoconductive devices can be characterized by their ability to form new bone tissue in non-osseous (soft-tissue) sites.

- Donah1 H, Kushitani S, Yasukawa E, Iwaki H, Hench LL, Wilson J, Tsai E, Sughran T: Particulate Bioglass Compared with Hydroxyapatite as a Bone Graft Substitute. *Clin Orthop* 334:316-325, 1997.
- Fujiyoshi Y, Hench LL, Corsiari H. Quantitative Rates of *In Vivo* Bone Generation for Bioglass and Hydroxyapatite Particles as Bone Graft Substitutes. *J Mater Sci Mat Med* 8:649-652, 1997.
- Vrounevelder WCA, Groot CG, de Groot K: Histological and Biochemical Evaluation of Osteoblasts Cultured on Bioactive Glass, Hydroxyapatite, Titanium Alloy and Stainless Steel. *J Biomed Res* 27:465-475, 1993.
- Yinos ID, Hukuman MW, Batten JJ, Buttery LD, Hench LL, Polak JM: Bioglass 45S5 Stimulates Osteoblast Turnover and Enhances Bone Formation *In Vitro*: Implications and Applications for Bone Tissue Engineering. *Calcif Tissue Int* 67:321-329, 2000.
- Yinos ID, Edgar AJ, Buttery LDK, Hench LL, Polak JM: Ionic Products of Bioactive Glass Dissolution Increase Proliferation of Human Osteoblasts and Induce Insulin-Like Growth Factor I mRNA Expression and Protein Synthesis. *Biochem Biophys Res Comm* 276:461-465, 2000.
- Boettch M, Cannas M: The Effect of Bioactive Glasses on Bone Marrow Stromal Cells Differentiation. *Biomaterials* 26(18):3873-3879, 2005.

#### HANDLING & USE






KINEX® implants should be combined with autogenous bone marrow aspirate at a ratio of 1:1 (KINEX®:BMA). KINEX® should not be used alone.

KINEX® should be implanted into bony defects according to the following technique. Prepare the wall of the defect with the contact of the KINEX® product, as needed. Mix or saturate the KINEX® product with autogenous bone marrow aspirate. Gently pack the site but avoid overfilling the bone void or compressing the treatment site. Remove excess material from the treatment site. Close the site using standard closure techniques and discard any unused KINEX® product.

**STORAGE**  
Do not expose to extreme heat. Store between 5°C (41°F) and 25°C (77°F).

**CONTACT INFORMATION**  
Globus Medical may be contacted at 1-866-GLOBUS1 (456-2871). A surgical technique manual may be obtained by contacting Globus Medical.

**STERILIZATION**  
KINEX® Bioactive is provided sterile in an unopened and undamaged package. This product must be used on or before the expiration date that is provided on the package label. This device is for single patient use and should never be reused.

SYMBOL TRANSLATION			
	CATALOGUE NUMBER		STERILIZED BY ETHYLENE OXIDE
	LOT NUMBER		STERILE BARRIER SYSTEM
	USE BY (YYYY-MM-DD)		CONSULT INSTRUCTIONS FOR USE
	CAUTION		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	SINGLE USE ONLY		AUTHORISED REPRESENTATIVE IN SWITZERLAND
	DO NOT RESTERILIZE		MANUFACTURER
	DO NOT USE IF PACKAGE IS DAMAGED		CONTAINS BIO ANIMAL
	PRODUCT QUANTITY		TEMPERATURE RANGE
	MEDICAL DEVICE		PRESCRIPTION USE ONLY

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