GLOBUS
MEDICAL

GLOBUS MEDICAL, INC.
1000 General Armitstead Avenue
Auburn, PA 19403
USA
Customer Service: 1-866-GLOBUS3 (101)
Phone: 1-866-466-2871
Fax: 1-866-466-2873

IMPORTANT INFORMATION ON THE HAVEN™ LAMINOPLASTY FIXATION SYSTEM
INFORMACIÓN IMPORTANTE Sobre el SISTEMA DE FIJACIÓN PARA LAMINOPLASTICA HAVEN™
INFORMAZIONI IMPORTANTI SUL SISTEMA DI FISSAZIONE PER LAMINOPLASTICA HAVEN™
INFORMACIONES IMPORTANTES SOBRE EL SISTEMA DE FIJACIÓN PARA LAMINOPLASTICA HAVEN™
INFORMAÇÕES IMPORTANTES SOBRE O SISTEMA DE FIXAÇÃO PARA LAMINOPLASTICA HAVEN™
BELANGRIKE INFORMATIE OVER HAVEN™ LAMINOPLASTIC FIXATIESSTEL
VÍTIGE INFORMATION OM HAVEN™ FERINGSSTYLLEN TIL LAMINOPLASTIK
VÍTIGE INFORMATION OM HAVEN™ FERINGSSTYLLEN TIL LAMINOPLASTIK
TÄVLIG INFORMATION OM HAVEN™ FERINGSSTYLLEN TIL LAMINOPLASTIK
DARWININGE INFORMATION PÅ HAVEN™ STYLLESTYLLEN TIL LAMINOPLASTIK HAVEN™

[ECB#] AJW Technology Consulting GmbH

Kongialallee 106, 40215 Düsseldorf, Germany

AUSTRALIA SPONSOR: GLOBUS MEDICAL AUSTRALIA PTY LIMITED, Unit

9/5-7 Inglewood Place Beaumaris Hills NSW 2153 Australia

0297

FOR OUTSIDE THE US ONLY

ENGLISH IMPORTANT INFORMATION
ON THE HAVEN™ LAMINOPLASTY FIXATION SYSTEM

DESCRIPTION
The HAVEN™ Laminoplasty System consists of spinal fixation plates for use in laminoplasty. These implants are manufactured from titanium or titanium alloy, as specified in ASTM F67, F136, F1295 and F1472. The plates may be used with previously cleared CANOPY®, ELLIPSE®, PROTEX CT®, QUARTEX®, RELIEVE® screws and CANOPY® Spacers.

INDICATIONS
The HAVEN™ Laminoplasty Fixation System is intended for use in the lower cervical and upper thoracic spine (C3-T3) in laminoplasty procedures. The HAVEN™ Laminoplasty Fixation System is used to hold bone graft material in place in order to prevent the graft from expulsing or impinging the spinal cord.

WARNINGS
One potential risk identified with this system is death. Other potential risks which may require additional surgery, include:

- device component fracture;
- loss of fixation;
- non-union;
- fracture of the vertebrae;
- neurologic injury; and
- vascular or venous injury.

Certain degenerative diseases or underlying physiological conditions such as diabetes, rheumatoid arthritis, or osteoporosis may affect the healing process, thereby increasing the risk of implant breakage or spinal fracture.

Components of this system should not be used with components of any other manufacturer.

The components of this system are manufactured from titanium or titanium alloy. Mixing of stainless steel implant components with different materials is not recommended for metallurgical, mechanical and functional reasons.

PRECAUTIONS

The implantation of laminoplasty fixation devices should be performed only by experienced spinal surgeons with specific training in the use of this system because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Preoperative planning and patient anatomy should be considered when selecting implant size.

Adequately sterilize the patient. Mental or physical impairment which compromises or prevents a patient's ability to comply with necessary limitations or precautions may place that patient at a particular risk during postoperative rehabilitation.

Surgical implants are SINGLE USE ONLY and must never be reused. An explanted implant must never be reimplanted. Even though the device appears undamaged, it may have small defects and internal stress patterns which could lead to breakage.

The HAVEN™ polyaxial screw hole plate is only to be used with ELLIPSE® polyaxial screws or QUARTEX® polyaxial screws. The HAVEN™ large polyaxial screw hole plate is only to be used with QUARTEX® polyaxial screws or QUARTEX® polyaxial plates. When using the HAVEN™ Laminoplasty Fixation System, the surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc., which may impact on the performance of this system.

CONTRAINDICATIONS

1. Use of the HAVEN™ Laminoplasty Fixation System is contraindicated when there is active systemic infection, infection localized to the site of the proposed implantation, or when the patient has demonstrated allergy or foreign body sensitivity to any of the implant materials.

2. Severe osteoporosis may prevent adequate fixation and thus preclude the use of this or any other orthopedic implant.

3. Components of this system are placed in close proximity to bone and implants, such as severe obesity or degenerative diseases, are relative contraindications. The decision whether to use these devices in such patients must be made by the physician taking into account the risks versus the benefits to the patient.

4. Use of these implants is relatively contraindicated in patients whose activity, mental capacity, mental illness, alcoholism, drug abuse, occupation, or lifestyle may interfere with their ability to follow postoperative restrictions and who may place undue stress on the implant during bony healing and may be at a higher risk of implant failure.

PACKAGING

The implants and instruments may be supplied pre-packaged and sterile, using gamma irradiation. The integrity of the sterile packaging should be checked to ensure that the sterility of the contents is not compromised. Packaging should be carefully checked for completeness and all components should be carefully checked to ensure that there is no damage prior to use. Damaged packages or products should not be used, and should be returned to Globus Medical. During surgery, after the correct size has been determined, remove the product from its sterile packaging.

The instrument sets are provided nonsterile and are sterile sterilized prior to use, as described in the STERILIZATION section below. Reusing or exposing to soil, instruments must be cleaned, as described in the CLEANING section below.

HANDLING

All instruments and implants should be treated with care. Improper use or handling may lead to damage and/or possible malfunction. Products should be checked to ensure that they are in working order prior to surgery. All products should be inspected prior to use to ensure that there is no nonacceptable deterioration such as corrosion, discoloration, pitting, cracked seals, etc. Non-working or damaged instruments should not be used, and should be returned to Globus Medical.

CLEANING

All instruments that can be disassembled must be disassembled for cleaning. All handles must be detached. Instruments may be cleaned following sterilization. The instruments should be cleaned using neutral cleaners before sterilization and introduction into a sterile surgical field or (if applicable) return of the product to Globus Medical.

Cleaning and disinfecting of instruments can be performed with alcohol-free solvents at higher temperatures. Cleaning and decontamination must include the use of neutral cleaners followed by a deionized water rinse. Note: certain cleaning solutions, such as those containing formalin, glutaraldehyde bleach and/or other alkaline cleaners may damage some devices, particularly instruments; these solutions should not be used.

The following cleaning methods should be observed when cleaning instruments after use or exposure to soil, and prior to sterilization:

1. Immediately following use, ensure that the instruments are wiped down to remove all visible soil and kept from drying by submerging or covering with a wet towel.

2. Disassemble all instruments that can be disassembled.

3. Rinse the instruments under running tap water to remove all visible soil. Flush the lumens a minimum of 3 times, until the lumens flush clean.

4. Prepare Enzol® (or a similar enzymatic detergent) per manufacturer's recommendations.

5. Immerse the instruments in the detergent and allow them to soak for a minimum of 2 minutes.

6. Use a soft bristled brush to thoroughly clean the instruments. Use a pipe cleaner for any lumens. Pay close attention to hard to reach areas.

7. Using a dry syringe, draw up the enzymatic detergent solution. Flush any lumens and hard to reach areas until no soil is seen exiting the tube.

8. Remove the instruments from the detergent and rinse them in running warm tap water.

9. Prepare Enzol® (or a similar enzymatic detergent) per manufacturer's recommendations.

10. Completely immerse the instruments in the ultrasonic cleaner and ensure detergent is in lumen by flushing the lumens.

11. Rinse the instruments from the detergent and rinse them in running deionized water or reverse osmosis water for a minimum of 2 minutes.

12. Dry instruments using a clean soft cloth and filtered pressurized air.

13. Visually inspect each instrument for visible soil. If visible soil is present, then repeat cleaning process starting with Step 3.

CONTACT INFORMATION

Globus Medical may be contacted at 1-866-GLOBUS1 (456-2871). A surgical technique manual may be obtained by contacting Globus Medical.

STERILIZATION

These implants are available sterile or nonsterile. Instruments are available nonsterile.

Stainless implants are sterilized by gamma radiation, validated to ensure a Sterility Assurance Level (SAL) of 10⁻⁶. Sterile products are packaged in a heat sealed, Tyvek pouch. The expiration date is provided in the package label. These products are considered sterile unless the packaging has been opened or damaged.

Nonsterile implants and instruments have been validated to ensure an SAL of 10⁻⁶. Use of a sterile is recommended, per the Association for the Advancement of Medical Instrumentation (AAMI) ST79, Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities. It is the end user's responsibility to use only sterilizers and accessories (such as sterilization wrap, sterilization pouches, chemical indicators, biological indicators, and sterilization cassette) that are designed for the selected sterilization cycle conditions (time and temperature). Sterile implants meet pyrogen limit specifications.

When using a steam sterilization container, the following must be taken into consideration for proper sterilization of Globus devices and loaded graphic cases:

• Recommended sterilization parameters are listed in the table below.

• Only rigid sterilization containers for use with pre-vacuum steam sterilization may be used.

• When selecting a rigid sterilization container, it must have a minimum filter area of 176 in² total, or a minimum of four (4) 7.5 in diameter filters.

• No more than one (1) loaded graphic case or its contents can be placed directly into a rigid sterilization container.

• Stand-alone modules/racks or single devices must be placed, without stacking, in a container basket to ensure optimal ventilation.

• The rigid sterilization container manufacturer's instructions for use are to be followed; if questions arise, contact the manufacturer of the specific container for guidance.

• Refer to AAMI ST79 for additional information concerning the use of rigid sterilization containers.

For implants and instruments provided NONSTERILE, sterilization is recommended (wrapped or contained) as follows:

Method	Cycle Type	Temperature	Exposure Time	Drying Time
Steam	Pre-vacuum	132°C (270°F)	4 minutes	30 minutes

These parameters are evaluated to sterilize only this device. If other products are added to the sterilizer, the recommended parameters are not valid and new cycle parameters must be established by the user. The sterilizer must be properly installed, maintained, and calibrated. Ongoing testing must be performed to confirm inactivation of all forms of viable microorganisms.

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